

VET-2 SONOTRON

"For results you can see in Training and count on at the Finish Line"

The Modality Chosen
By Trainers for
keeping their Horses
in the Winners Circle

VET-2 Sonotron
Therapy can be used
in conjunction with
your daily training
regiment and does
not alter
conventional
veterinary
treatments.

Great Adjunctive
Treatment for
Pre-Race
Conditioning

Pain, Inflammation
Tissue Healing
Lameness

*"Keeping You in Front of the
Line"*



VET-2 SONOTRON THERAPY

Therapeutic Technology for Horses

VET-2 SONOTRON THERAPY

International Distributor:
Panmedic Sdn Bhd,
A-3-12, 8 Avenue, Jalan Sungai
Jernih 8/1, 46050 Petaling Jaya,
Selangor, Malaysia.
Tel: (603) 79543622
Fax: (603) 79541622
Websites: www.panmedic.com.my
www.sonotron.com.my
Email: vet@sonotron.com.my

Rapid
Noninvasive
Without Drugs
Without Side
Effect
Excellent
Animal
Compliance



VET-2 SONOTRON

VET-2 SONOTRON THERAPY uses the most advanced technology for treating joint problems and inflammation in horses. This technology was proven to reduce lameness and increase range of motion for horses in controlled laboratory studies and on the track by leading trainers within the industry. The therapy has no harmful side effects, is non-invasive and does not cause pain to the horses. It is the safe, efficient way to reduce lameness and improve the condition of horses with joint problems .

- Proven on the track by Top Trainers
- Boost Performance Prior to Race Day
- Safely used along with Conventional Treatment
- Can be easily used by grooms in their regular training schedule

INDICATIONS FOR USE

Pain, Inflammation
Tissue Healing
Lameness

FOR TREATING

Jacks, Stifles, Hocks, Knees,
Feet, Tendons, Body

What Trainers Have Reported...

"Results after the First Treatment"

"Achieved Results while Dramatically Cutting Costs"

"Ease of use" - grooms can be quickly trained to administer treatments



Clinically Tested

Extensive double-blind studies were conducted on ponies with induced arthritis at the University of Wisconsin Veterinary School. The results of the studies showed that the technology used by VET-2 SONOTRON reduced the level of joint problems in the horses and increased their range of motion. Swelling was reduced and gait analysis revealed a marked improvement in range of motion. The results of this study were published in the well-respected Canadian Journal of Veterinary Research in early 1991.

Results at the Track

Horses treated with VET-2 SONOTRON will perform better. Two outstanding examples of horses that have responded well to the treatment have been Jet Lair and Fine Cuisine, both out of Freehold, New Jersey. Their trainer has confirmed that these horses have significantly improved since starting treatment. This has resulted in top finishes in a number of upper level stakes races. The most recent races won by Fine Cuisine are The All Pennsylvania Sire Stakes and The New Faces Race. Jet Lair has posted very impressive finishes, including 4/10/98 finishing 1st, 4/17/98 - 2nd, 5/16/98 - 1st, 5/22/98 - 3rd, and 5/29/98 - 2nd .

Wajastar, a 3-year-old was diagnosed as having degenerative joint disease syndrome or osteoarthritis in the left rear hock with some secondary sensitivity in the back. This was determined by experts at the University of Pennsylvania School of Veterinary Medicine, New Bolton Center on June 3, 1991. The Center recommended 3 to 4 months of rest before resuming racing. VET-2 SONOTRON THERAPY was started on Wajastar on July 15, 1991. On 7/31/91 Wajastar finished 2nd. On 8/8/91 he finished 4th. On 8/17/91 he finished 2nd losing by a head but running one of the fastest times of his career. His trainer

reported that not only did he race well but all signs of the lameness he was showing prior to treatment disappeared.

These are results you can expect from VET-2 SONOTRON THERAPY.

How Does It Work?

VET-2 SONOTRON THERAPY employs Pulsed Radio Frequency Therapy or PRFT. The device produces a safe corona discharge beam of PRFT, which is administered to the affected joint for a certain period of time over several treatment visits and can be used on the day of the race with positive results. PRFT has no harmful side effects but has been proven to reduce lameness and pain in inflamed joints and soft tissues.

The Treatment Regimen

Once a horse is determined to have a pain related problem, VET-2 SONOTRON THERAPY should be commenced. The standard treatment regimen recommended entails 3 to 4 treatment visits over a 2 to 4 week period with treatment visits a minimum of 2 days apart. At each visit the horse will have a series of Treatment Units (TUs) applied to the area which can be performed in the stall or other convenient location. After treatment, the horse is walked for approximately 3 minutes. The entire treatment visit takes approximately 15- 20 minutes. Other than this simple procedure, there is no special handling required and no change in training necessary. Grooms can easily incorporate VET-2 SONOTRON THERAPY into their daily training routines.

Results should be observed after the first or second visit and maximum effect realized after the third or fourth. After the initial treatment regimen is completed it is advisable to treat again prior to each race for maintenance.

Anyone interested in the VET-2 Sonotron please contact Panmedic Sdn Bhd, A-3-12, 8 Avenue, Jalan Sungai Jernih 8/1, 46050, Petaling Jaya, Selangor, Malaysia. Tel: (603) 79543622, Fax: (603) 79541622, Email: vet@sonotron.com.my